

Gel documentation system



The Paya- gel documentation system allows users to perform image capture functions with user-friendly panel. The simple design is excellent for multiple user laboratories. It can easily connect to computer by though USB cable to transfer images directly to your computer. The reliable and safe Paya gel documentation is an excellent choice for molecular biology laboratories.



Technical specifications

UV TRANSILLUMINATOR	256 nm & 365nm
FILTER	580AF12062mm Ethidium Bromide Filter
Resolution	18 Megapixels
Lens	55-250mm
Data type	JPEG
UV Transilluminator size	20*20 cm
Electrical Requirement	220 V AC - 50HZ
WEIGHT	10 KG
DARKROOM DIMENSION (W*H*D)	34*67*44 CM